


## Wie Sie Ihre Kaleido/Diabeloop DBLG1 Daten hochladen

# Kaleido.




Für den Datenexport

1. Öffnen Sie die Seite „[Anmelden | YourLoops](#)“ und melden Sie sich mit Ihren persönlichen Zugangsdaten an.
2. Klicken Sie oben rechts in der Menüleiste auf den Button «Download-Berichte».



3. Wählen Sie als Zeitraum den Zeitraum Ihres letzten Menstruationszyklus aus und exportieren Sie die Dateien im CSV-Format.
4. Starten Sie den Export und speichern Sie die heruntergeladene .zip-Datei auf Ihrem Computer.
5. Entpacken Sie die .zip-Datei und laden Sie die enthaltenen Dateien über den dafür vorgesehenen sicheren Link hoch.



The screenshot shows the 'Download-Bericht' form. It has two main sections: 'Festen Zeitraum wählen' (Choose fixed period) and 'oder benutzerdefinierter Zeitraum' (or user-defined period). Under 'Festen Zeitraum wählen', there are buttons for '1 Woche', '2 Wochen', '4 Wochen', and '3 Monate'. Under 'oder benutzerdefinierter Zeitraum', there is a date range selector showing '12. Sept. 2025 → 18. Sept. 2025'. Below this is a calendar view for August 2025 and September 2025. At the bottom, there is a section 'Ausgabeformat wählen' (Choose output format) with radio buttons for 'PDF' (selected) and 'CSV'. At the very bottom, there are two buttons: 'Abbrechen' (Cancel) and 'Herunterladen' (Download).